



PHYSIO
RADIANCE

Visage⁺



Physio Radiance Visage+ Facial Device

THE EVOLUTION OF BEAUTY



GROWING YOUNGER
GLOWING MORE RADIANT
THAN EVER

#StartGlowingNow

INTRODUCTION



A revolutionary facial Treatment device specially designed to help you look younger and more radiant.

This device is a breakthrough **anti-ageing** treatment that brings together **4 types of technology** to make your **skin glow** like never before

INTRODUCTION



Integrating the latest in **skincare technology** and trends, this device comes in a small, sleek, easy-to-carry package, with simple functions that **rejuvenate your skin at your convenience.**

For optimal results, recommend to use the device **2-3 times a week.**

KEY BENEFITS

- Delay the signs of ageing by reducing fine lines and wrinkles.
- Lighten age spots.
- Enhance skin hydration levels.
- Improve skin texture.
- Tone facial muscles.



FUNCTION INSIGHT : 4 in 1 system



PHYSIO RADIANCE VISAGE+
provides the most advanced, safe and
non-invasive treatment by harmonising
4 types of technology.

1st : ELECTRICAL MUSCLE STIMULATION



1



ELECTRICAL MUSCLE STIMULATION (EMS)

EMS helps tighten your skin, promotes facial elasticity, as well as gives you a facelift. It also increases muscle strength by improving skin firmness and tone, ultimately reducing facial fat.

2nd : RADIO FREQUENCY



RADIO FREQUENCY
(RF)

RF helps smoothen fine lines by stimulating facial tissue to produce new collagen and elastin.



3rd : CHROMO THERAPY



CHROMO THERAPY

Light Therapy is great for calming irritation and brightening the skin. It's also a brilliant approach to trigger your body's natural cell processes, to accelerate rejuvenation and repair the skin. It is recommended for skin rejuvenation, treating sun damage, fighting acne-causing bacteria, inflammatory conditions, scarring, plumping the skin, and reducing wrinkles.

3rd : CHROMO THERAPY




Red Light

- Improves firmness
- Stimulates the production of collagen and elastin
- Increases blood circulation



Green Light

- Reduces skin pigmentation
- Brightens skin
- Reduces fine lines



Yellow Light

- Smoothens skin
- Reduces redness



Blue Light

- Improves acne and prevent breakouts
- Calms skin by reducing excessive sebum

4th : GALVANIC TECHNOLOGY



GALVANIC TECHNOLOGY


This technology helps deep cleansing process while providing intense hydration. It increases the effectiveness of Physio Radiance in moisturising your skin and also helps deliver the product's ingredients deep into your skin.



GETTING TO KNOW YOUR VISAGE+

Mode Button (M) :
To turn on/off
and select Mode

LED Button (L) :
To select LED Lights
(Red/Green/Yellow/Blue)



Mode & Battery
Indicator

- C : Cleansing
- M : Moisturising
- R : Radiance-inducing
- L : Lifting

Head USB Charging
Port


YOUR VISAGE+ KIT



The following items are provided in Physio Radiance Visage+ kit

1. Physio Radiance Visage+ Device
2. USB Cable
3. Physio Radiance Visage+ Travel Pouch
4. User Manual

QNET



For more information please **CLICK** on **ICON** above